

URIT

A better choice for portable & reliable
urinalysis with multifunction

URIT-31 Urine Analyzer



Global Diagnostics Supplier

URIT-31

Urine Analyzer

Specifications

Parameters: pH, Nitrite, Specific Gravity, Blood, Glucose, Bilirubin, Urobilinogen, Ketone, Leucocytes, Protein, Ascorbic Acid, Creatinine, Microalbumin, Calcium and urine color

Measurement Principle: dual wavelength reflectance method

Measuring Wavelength: 4 LED wavelengths (470, 550, 620, 720nm)

Sample preparation: manual dipping

Reaction time: 60 seconds per test

Data memory: up to 1000 samples

Display: LCD colored touch-screen

Strips: URIT 10G/11G/14G urine reagent strip.

Interface: Micro USB, Bluetooth

Work condition: environmental temperature 5–40°C, RH:20%~80%

Power supply: power adapter 9V DC, 1.5A

Dimensions: 120(L) × 83(W) × 31(H)mm (analyzer)
110(L) × 80(W) × 38(H)mm (printer)

Weight: 180g (analyzer), 180g (printer)



Features

- Point-of-Care urinalysis with palm-sized feather weight unit
- Bluetooth technology for wireless printing and data transfer
- LCD colored touch-screen
- Single button easy operation
- Recharge, transfer data and upgrade, all with Micro USB cable
- Compatible with URIT 10G, 11G and 14G urine reagent strips



URIT URINE CONTROL



URIT 10G、11G、14G urine reagent strip

14G with Micro albumin, Creatinine and Calcium.

URIT Medical Electronic Co., Ltd.

Head Office

Add.: No. D-07 Information Industry District,
High-tech Zone, Guilin, Guangxi 541004,
P. R. China

Tel: +86-773-2288586 Fax: +86-773-2288560

Shenzhen Office

Add.: Room E.J, 24th Floor, New Baohui Building, Nanshan Road,
Nanshan District, Shenzhen, P.R. China

Tel: +86-755-26050853 Fax: +86-755-86036150

E-mail: export@uritest.com service@uritest.com Http://www.urit.com



ISO9001
ISO13485