

# Cardiac Monitoring

## NIDEK ECG Machines

ECG-701 | ECG-703 | ECG-703A | ECG-712

Portable, Precise, User-friendly

Nidek India's Electrocardiograph machines are one of the basic diagnostic tools for cardiologists and general physicians.

They can be used to perform routine diagnostic checks. The machines are simple to use and precise, thus, making the diagnosis easy to perform and analyze.

Nidek's ECG machines are available in several models suited to the physicians and patient's needs and requirements.

Product descriptions appear in the following pages.



## ECG-701

### Single channel ECG machine



#### Key Features

- Single channel ECG machine with interpretation
- Colour display
- Manual / Auto / RR / Store Modes (optional)
- Built in rechargeable battery
- Adapt to 110-230V, 50 / 60Hz AC power supply
- Built in RS232 interface (ECGNET software optional)

#### Specifications

Display	LCD 3.5"
Paper	50mm X 20m roll / fold paper
Dimension	285 X 265 X 72 Unit (mm)

## ECG-703

### Three channel ECG machine

#### Key Features

- 3 channel ECG machine with interpretation
- Colour display
- Manual / Auto / RR / Store Modes
- 80mm 3-ch recording with rhythm lead
- Built in rechargeable battery
- Adapt to 110-230V, 50 / 60Hz AC power supply
- Real-time ECG waveform, cases save / replay
- Built in USB / RS232 interface (ECGNET software optional)



#### Specifications

Display	LCD 3.5"
Paper	80mm X 20m roll / fold paper
Dimension	285 X 265 X 72 Unit (mm)

## ECG-703A

### Three channel ECG machine

#### Key Features

- 3 channel ECG machine with interpretation
- Colour display
- Manual / Auto / RR / Store Modes
- 80mm 3-ch recording with rhythm lead
- Built in rechargeable battery
- Adapt to 110-230V, 50/60Hz AC power supply
- Real-time ECG waveform, cases save/replay
- Built in USB / RS232 interface (ECGNET software optional)



#### Specifications

Display	LCD 5.7"
Paper	80mm X 20m roll / fold paper
Dimension	285 X 265 X 72 Unit (mm)

## ECG-712

### Twelve channel ECG machine

#### Key Features

- Portable 12 channel ECG machine
- Colour display
- Manual / Auto / RR / Store Modes
- 7 inch LCD
- Alphabetic keyboard for patient information
- Built in rechargeable battery
- Adapt to 110-230V, 50 / 60Hz AC power supply
- Real-time ECG waveform, cases save / replay
- Built in USB / RS232 interface (ECGNET software optional)



#### Specifications

Display	LCD 7.0"
Paper	210mm X 216mm / 20m roll / Z Fold paper
Dimension	320 X 270 X 66 Unit (mm)

## Standard Specifications

Lead	Standard 12 leads
Acquisition	Simultaneous 12 leads
Resolution	12bit / 1000Hz
Operation mode	Manual / Auto / RR / Store Modes
Filter	AC, EMG Filter
Drift filter	Anti-Drift System: 0.15 Hz (Adaptive)
CMRR	$\geq 120\text{dB}$ , with AC Filter
Input circuit	Floating, Protection against defibrillator effect
Input circuit current	$\leq 0.05 \mu\text{A}$
Input Impedance	$> 50 \text{M}\Omega$
Patient current leakage	$< 10 \mu\text{A}$
Calibrating voltage	$1\text{mV} \pm 2\%$
Voltage tolerance	$\pm 500\text{mV}$
Time constant	$> 3.2\text{s}$
Frequency response	$0.05\text{Hz} \sim 160\text{Hz} (-3\text{db})$
Noise level	$\leq 15 \mu\text{V}_{\text{p-p}}$
Threshold	$\leq 20 \mu\text{V}$
Paper speed	6.25, 12.5, 25, 50mm / st ( $\pm 3\%$ )
Sensitivity	Auto, 2.5, 5, 10, 20, 40mm/mV
Recorder	High resolution thermal printer
Safety standard	GB9706.1-2007, IEC Class I, Type CF

\* Specifications subject to change without notice.



### Nidek Medical India Private Limited

**Head Office:** 238 B, A.J.C. Bose Road, Kolkata 700 020  
Tel: (033) 4017 7000, Mob: 93395 66444

**Delhi:** 322, 3rd floor, Ansal Chamber II, 6, Bhikaji Cama Place, New Delhi 110 066  
Tel: (011) 2616 5742, Mob: 93501 83222

**Chennai:** Plot No. 1966, J Block, 3rd Street, 13th Main Road, Anna Nagar (W), Chennai 600 040  
Tel: (044) 4266 3662, Mob: 93808 06777, 93459 70044

**Mumbai:** C-216, Kailash Ind. Compl., Veer Savarkar Marg, Hiranadani Link Road, Vikhroli (W), Mumbai 400 079  
Tel: (022) 4215 5441, Mob: 73039 11666

**Hyderabad:** House Bearing No. 13/1, Shiva Arun Colony, West Marredpally, Secunderabad 500 026  
Tel: (040) 4015 0243, Mob: 93468 32266

**Lucknow:** 1/44, Vijayant Khand Gomti Nagar, Lucknow 226 010  
Mob: 93361 20941

**Chandigarh:** 2nd Floor, SCO# 355, Sector 32-D, Chandigarh 160 031  
Mob: 93562 78899

E-mail: [sales@nidekindia.com](mailto:sales@nidekindia.com), Website: [www.nidekindia.com](http://www.nidekindia.com)